

IPD Implant Level Angulation Guide Chart

Brand:		Nobel Biocare (SERIE A)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Branemark System	Ø3,5	Available	30°	Available	30°
	Ø4,1	Available	30°	Available	30°
	Ø5,1	Available	30°	Available	30°
Multiunit	Ø4,8	Available	30°	Available	30°
Active	Ø3,0	Available	0°	Available	0°
	Ø3,5	Available	30°	Available	30°
	Ø4,3	Available	30°	Available	30°
Replace select	Ø3,5	Available	20°	Available	20°
	Ø4,3	Available	25°	Available	25°
	Ø5,0	Available	25°	Available	25°

Brand:		Biomet 3I (SERIE B)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Externa Osseotite	Ø3,4	Available	30°	Available	30°
	Ø4,1	Available	30°	Available	30°
	Ø5,0	Available	30°	Available	30°
Certain	Ø3,4	Available	25°	Available	25°
	Ø4,1	Available	30°	Available	30°
	Ø5,0	Available	30°	Available	30°

Brand:		Klockner (SERIE C)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
SK2/NK2	Ø4,2	Available	30°	Available	30°
Essential Cone	Ø4,5	Available	30°	Available	30°
KL	Ø3,5	Available	30°	Available	30°
	Ø4,1	Available	30°	Available	30°
	Ø5,1	Available	30°	Available	30°
Vega	Ø3,5	Available	30°	Available	30°
	Ø4,1/Ø4,5	Available	30°	Available	30°

IPD Implant Level Angulation Guide Chart

Brand:		Straumann (SERIE D)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Tissue Level	Ø4,8	Available	30°	Available	30°
	Ø6,5	Available	30°	Available	30°
Bone Level	Ø3,3 NC	Available	25°	Available	25°
	Ø4,1 RC	Available	25°	Available	25°
BLX	RB/WB	Not available	0°	Not available	0°

Brand:		Astra (SERIE E)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Osseospeed	Ø3,0	Available	0°	Available	0°
	Ø3,5/Ø4,0	Available	25°	Available	25°
	Ø4,5/Ø5,0	Available	25°	Available	25°
Evolution	Ø3,6	Available	30°	Available	30°
	Ø4,2	Available	25°	Available	25°

Brand:		Zimmer (SERIE F)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
TSV	Ø3,5	Available	30°	Available	30°
	Ø4,5	Available	30°	Available	30°
	Ø5,7	Available	30°	Available	30°
Swissplus	Ø4,8	Available	30°	Available	30°
Eztetic	Ø3,1	Available	30°	Available	30°

Brand:		BTI (SERIE G)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Externa	Ø3,5	Available	30°	Available	30°
	Ø4,1	Available	30°	Available	30°
Interna Universal	Ø4,1	Available	30°	Available	30°
Multi-IM	Ø4,1	Available	30°	Available	30°

IPD Implant Level Angulation Guide Chart



Brand:		Microdent (SERIE H)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Universal	Ø4,1	Available	30°	Available	30°
	Ø5,0	Available	30°	Available	30°

Brand:		Dentsply (SERIE I)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Xive	Ø3,4	Available	30°	Available	30°
	Ø3,8	Available	30°	Available	30°
Ankylos	C/X	Not available	0°	Available	0°

Brand:		Camlog (SERIE J)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Camlog	Ø3,3	Not available	0°	Available	0°
	Ø3,8	Not available	0°	Available	0°
	Ø4,3	Not available	0°	Available	0°

Brand:		Biohorizons (SERIE L)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Tapped Internal	Ø3,0	Available	30°	Available	30°
	Ø3,5	Available	30°	Available	30°
	Ø4,5	Available	30°	Available	30°
	Ø5,7	Available	30°	Available	30°

IPD Implant Level Angulation Guide Chart



Brand:		Sweden & Martina (SERIE M)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Outlink	Ø3,3	Available	30°	Available	30°
	Ø4,1	Available	30°	Available	30°
	Ø5,0	Available	30°	Available	30°
Premium Kohno	Ø3,3	Available	30°	Available	30°
	Ø3,8	Available	30°	Available	30°
	Ø4,25	Available	30°	Available	30°

Brand:		Phibo (SERIE N)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
TSH	S2	Available	30°	Available	30°
	S3-S4	Available	30°	Available	30°
	S5	Available	30°	Available	30°

Brand:		Osstem Implants (SERIE O)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
TSIII	Ø3,5	Available	20°	Available	20°
	Ø4,0	Available	25°	Available	25°

IPD Implant Level Angulation Guide Chart



Brand:		Biotech Dental (SERIE Q)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Kontakt	RP		0°		0°

Brand:		Neodent (SERIE R)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
GM	GM		20°		20°
Helix	∅3,3		25°		25°
	∅4,1		30°		30°
	∅5,0		30°		30°
GM Abutment	RP		30°		
GM Micro	∅3,5		30°		30°

Brand:		Bego (SERIE S)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Semados	∅3,25/∅3,75		30°		30°
	∅4,1		30°		30°
	∅4,5		30°		30°

Brand:		Mis (SERIE T)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Seven	Narrow 3,3		30°		30°
	Standard		30°		30°
	Wide		30°		30°
C1/V3	Standard		30°		30°

IPD Implant Level Angulation Guide Chart

Brand:		Dio (SERIE U)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
UFII	NP	Available	30°	Available	30°
	RP	Available	20°	Available	20°

Brand:		Megagen (SERIE W)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Anyone	RP	Available	10°	Available	10°
Anyridge	Ø4,0	Not available	0°	Not available	0°

Brand:		Dentium (SERIE X)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
Implantium	RP	Available	10°	Available	10°

Brand:		Medentis (SERIE Y)			
System	Prosthetic Platform	Engaging		Non-Engaging	
		Availability	ASC max degree	Availability	ASC max degree
ICX	RP	Available	30°	Available	30°